

Cyanacrylate glue

Cyanacrylate glue is a rapid-hardening, solvent-free single component adhesive. This instant adhesive is used to bond different kinds of earmold materials.

Unit: Cyanacrylate glue 10 g

REF: 55200

EUH202! Cyanoacrylate. Danger. Bonds skin and eyelids in seconds. Keep out of reach of children! Contains: Ethyl-2-cyanoacrylate

Usage:

I. Close container properly after usage.

2. Avoid contact with skin. Before handling, protect hands with appropriate gloves.

Handling:

Apply Cyanacrylate glue always very sparingly. The thinner the bonded joint, the firmer the bond and the more quickly hardening takes place. The ultimate strength of the bond is achieved within 12 hours. The usage of Cyanacrylate glue requires handling by qualified personnel considering the relevant safety data sheet. Working with Cyanacrylate glue should be done under a vent. Please see safety data sheet of this product for disposal.

Handling instruction:

I. Preparation

Best results are achieved on clean surfaces. Remove dust and harder particles by using compressed air. It is advisable to clean the surface before-hand with *Surface Cleaner (REF 46400)* from greasy and oily remainders. Please consider the evaporation times of the used cleaner. The adhesive and the substrates to be bonded should be at a temperature between +18 °C and +25 °C before processing to optimize the adhesive strenght.

2. Gluing

Apply a thin film of Cyanacrylate glue. Depending on the material and the amount of glue used, press the parts to be glued together tightly. Due to the capillary action the adhesive penetrates the smallest gaps.

Safety regulations:



Warning!

H: 315-319-335: Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation. P: 102-261-302+352-305+351+338-312-501

Keep out of reach of children. Avoid breathing dust/fume/gas/mist/vapors/spray. IF ON SKIN: Wash with plenty of water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER/doctor if you feel unwell. Dispose of contents/container according to local and national requirements/regulations. Contains: Ethyl 2-cyanoacrylate.

Please follow the product information and the material safety data sheet





HEBA-OTOPLATIK Labortechnik u. Einrichtungs GmbH & Co. KG Hauptstr. 42 - D-63853 Mömlingen Tel: + 49 6022 681600 -Fax: + 49 6022 31663 info@heba.de -www.HEBA.de

Instruction for use

Applications:

This instant adhesive is used to bond different kinds of earmold materials.

Storage and durability:

The shelf life is, at optimum storage temperature of +2 °C to +10 °C in the closed original container, a maximum of 24 months. A higher storage temperature leads to a significantly lower durability. The storage temperature must not fall below +2 °C. Protect the product from direct sunlight or heat, frost and humidity.